

LIBLIA® H-FABP



H-FABP

LIBLIA® H-FABP is a reagent for quantification of H-FABP in serum or plasma based on latex-agglutination immunoassay method.

Heart-type fatty acid-binding protein (H-FABP) is a low-molecular-weight (15KDa) cytoplasmic protein involved in the intracellular uptake and buffering of free fatty acids in the myocardium¹⁾. When the myocardium is injured, H-FABP is easily released into the circulation and it can be detectable even if in super-acute phase within 2-4hrs after the onset of symptoms. Therefore, H-FABP is thought to be an excellent biomarker for acute myocardium damage that is Acute Coronary Syndrome (ACS), heart failure, etc.

LIBLIA® H-FABP is quick and convenient reagent based on latex-agglutination method, it's one of the most useful emergency laboratory test for detection of myocardial damage.



DS PHARMA BIOMEDICAL

DS Pharma Biomedical

DS Pharma Biomedical Co., Ltd.
33-94 Enoki, Suita
Osaka 564-0053, JAPAN
<http://www.dsp-bio.com>
TEL +81-6-6337-5941
FAX +81-6-6337-6020




Figure-1. ROC Curves of H-FABP, Myoglobin and CK-MB within 3hr after the onset of symptoms³⁾.

LIBLIA[®] H-FABP

Intended Use


For the quantitative measurement of H-FABP in serum or plasma.

Summary

H-FABP is sensitive and early marker for myocardial injury. Especially, H-FABP is famous as new marker of the Acute Coronary Syndrome (ACS) diagnosis, the potential of it as one of the most sensitive marker in super-acute phase within 2-4hr from onset²⁻⁶. H-FABP can be detectable in a sample from the patient who has small myocardial injury, therefore it is useful marker of acute heart failure patient's condition grasp, too⁷⁻⁹. Additionally, it has been reported that prognostic utility of H-FABP in patients with ACS in spite of with/without High Sensitive Troponins¹⁰⁻¹².

Principle

The LIBLIA H-FABP is an latex-agglutination turbidimetric immunoassay. Sample is added to a buffer solution and mixed with a suspension of mouse anti-human H-FABP monoclonal antibody that is bound to latex. H-FABP binds to the latex-bound antibody and agglutinates. The light scattering caused by the increase in particle size is used as a measure of H-FABP concentration. The amount of light scattering is proportional to the concentration of H-FABP in the sample.



Reagents

Composition

- R1-Reagent: Liquid, ready to use.
- R2-Reagent: Liquid, ready to use.

Storage and Expiry period

- Storage: Store in a cool place (2-10°C), protected from light. Avoid freezing.
- Expiry period: 1 years

Specimen Collection and Preparation

Serum, EDTA-plasma, and sodium or lithium heparinized-plasma are the recommended collection media. If not analyzed within 24hrs, they may be preserved at -20°C or below.

Assay Procedure

Figure-1 is a general example of the LIBLIA H-FABP assay procedure for an automated analyzer. All analyzer applications should be validated.




Figure-1

Performance

Precision

When two distinct samples are determined 5 times simultaneously, the coefficient of variation in their concentration should be less than 10 %.

Accuracy

Control samples should show a value within ±10% of its known concentration with this reagent.

Correlative Test (LIBLIA vs ELISA)

Comparative performance studies were conducted using LIBLIA H-FABP and Markit[™]-M H-FABP (ELISA). 104 Serum samples, with H-FABP concentrations between 2.3 and 106.8 ng/mL were tested (Figure-2).




Figure-2

References

- 1) Schaap FG., et al.:Mol.Cell Biochem., 180:43, 1998
- 2) Tanaka T., et al.:Clin.Biochem., 24:195, 1991
- 3) Okamoto F., et al.:Clin.Chem.Lab.Med., 38:231, 2000
- 4) Seino Y., et al.:Am.J.Med., 115:185, 2003
- 5) Nakata N., et al.:Cardiology, 99:96, 2003
- 6) McCann CJ., et al.: Eur.Heart J., 29(23):2843, 2008
- 7) Setsuta K., et al.:Am.J.Med., 113:717, 2002
- 8) Sato Y., et al.:Heart., 90:1110, 2004
- 9) Setsuta K., et al.:Circ.J., 72:569, 2008
- 10) Kilcullen N., et al.:J.Am.Coll.Cardiol., 50:2061, 2007
- 11) O'Donoghue M., et al.:Circulation, 114:550, 2006
- 12) Suzuki M., et al.: Int.Heart J.46:601, 2005